

1,000x Ca⁺⁺/Mg⁺⁺ Stock Solution for HBSS

(Greg A. Perry, Ph.D.)

(100 ml at 1,000x)

Equipment:

Beaker
Stir plate & stir bar
Sterile 100ml bottle
0.2 µm Vacuum Filter

Reagents:

100 ml	Distilled Water
14 gm	Calcium Chloride, anhydrous (CaCl ₂)
10 gm	Magnesium Chloride, hexahydrate (MgCl ₂ • 6H ₂ O)
10 gm	Magnesium Sulfate, heptahydrate (MgSO ₄ • 7H ₂ O)

Method:

1. Dissolve all solids in distilled water.
2. Sterile filter into glass container.

Use: Add 0.1 ml of stock for each 100ml of HBSS