

10x FACSLyse Solution

(Greg A. Perry, Ph.D.)

Reagents:

| | |
|-------|-------------------------|
| 2.5 g | Sodium Citrate |
| 27 ml | Formaldehyde (37%) |
| 30 ml | Diethylene Glycol (DEG) |
| 43 ml | Distilled Water |

Method:

- 1) Dissolve sodium citrate in distilled water.
- 2) Add formaldehyde and mix.
- 3) Add diethylene glycol. DEG may be difficult to get into solution, so mix well.

To Use:

- 1) Make a 1x working solution by diluting 10x FACSLyse solution with distilled water.
 - a. For example, add 1ml 10x stock to 9 ml distilled water.
- 2) Add 2ml of 1x working solution to 100µl of whole blood.
- 3) Vortex gently, and incubate 10 minutes at room temperature in the dark.
- 4) Centrifuge and wash twice with PBS.

Notes:

- B/D adds 370mg of Tetrasodium EDTA to the 10x stock, however it does not seem to be required.
- If using in conjunction with staining for flow cytometry, cells should be stained and washed prior to lysing. The lysing step should be done last, just before resuspending cells in FACSFix.