

Using the Coulter Counter Model Z1
(Detailed Instructions)

Startup (to be run at the beginning of each day)

- 1) Turn on the power switch.
 - The S1 setup screen will be displayed on the data terminal LCD.
 - 2) Check to make sure the upper and lower size limits and the count mode are appropriate. Generally:
 - $T_u = 21 \mu\text{m}$
 - $T_l = 6 \mu\text{m}$
 - Count Mode = Above T_l
 - 3) Press "Set-Up"
 - The S2 setup screen will be displayed on the data terminal LCD.
 - 4) Check to make sure that:
 - Select Aperture = 100 μm C
 - Metered Volume = 0.5
 - 5) Remove the vial containing detergent from the aperture, replace it with a vial containing 20 ml of fresh isotonic saline, and move it into position.
 - 6) Press the "Function" key.
 - The F1 function screen will be displayed on the data terminal LCD.
 - 7) Using the cursor keys (\leftarrow & \rightarrow) move to "Prime Aperture"
 - 8) Press the "Start / Stop" key to enter the priming routine.
 - 9) Press the "Start / Stop" key again to begin priming the aperture.
 - The aperture will be primed with saline. This takes about 5 minutes.
 - 10) Press "Output" to go to the output A1 page. The system is now ready to use.
-

Running Samples

- 1) Press the "Output" key until the Output A1 menu is displayed.
- 2) Under "Result Type" choose either:
 - Count ----- output displays the # of events counted.
 - Concentration ----- output displays the # of cells/ml, based on the dilution factor entered.
- 3) If you choose concentration, then enter the appropriate dilution factor using the number keys in the "Dilution Factor" field.

Count the background ...

- 4) Place a vial containing 10ml of fresh isotonic saline on the stage and move it into position.
- 5) Press the "Start / Stop" key to begin counting the background.
 - Background counts should be under 50 to continue.
 - If not, rinse the aperture and repeat steps 4 and 5.

Count your cell sample ...

- 6) Place a vial containing 10ml of isotonic saline and your cells on the stage and move it into position.
- 7) Press the "Start / Stop" key to begin counting.
- 8) Record the counts or cells/ml displayed on the LCD
 - Generally, each sample should be counted twice and a mean value taken to insure accuracy.
- 9) For your next sample, place the new sample on the stage and repeat steps 6-8.

(continued on next page)

Rinse the aperture ...

- 10) Place a vial containing 10ml of fresh isotonic saline on the stage.
 - 11) Rinse the aperture by moving the stage up & down over the aperture several times.
 - 12) Place a vial containing 10ml of fresh isotonic saline on the stage and move it into position.
 - 13) Press the "Start / Stop" key to begin counting the background.
 - Background counts should be under 50 to continue.
 - If not, rinse the aperture and repeat steps 10-13.
-

Shutdown *(to be run at the end of each day)*

- 1) Place a vial containing detergent (*Coulter Clenz*) on the stage and move it into position.
- 2) Press the "Function" key.
 - The F1 function screen will be displayed on the data terminal LCD.
- 3) Using the cursor keys (← & →) move to "Prime Aperture"
- 4) Press the "Start / Stop" key to enter the priming routine.
- 5) Press the "Start / Stop" key again to begin priming the aperture.
 - The aperture will be primed with detergent. This takes about 5 minutes.
- 6) Turn off the power switch.