

Immunocytochemistry Protocol

Greg A. Perry, Ph.D.

Equipment:

Frozen Tissue
Microtome

Gelatin coated glass slides
Coplin Jars

Reagents:

DPBS with 5% FCS and 0.05% Tween-20
Tris buffer
Antibodies
Sodium bicarbonate buffer
Diaminobenzadine (DAB) working solution
Osmium Tetroxide (OsO₄) working solution

1% Gluteraldehyde in DPBS
95% ETOH
100% ETOH
Xylene (or HistoClear)
Permount

Method:

- 1) Cut 6 µm sections of tissue; mount on gelatin coated slides.
- 2) Immerse 1-2 seconds in cold acetone (on ice).
- 3) Air dry slide 30-60 seconds.
- 4) Immerse 15-60 minutes in DPBS + 5% FCS + 0.05% Tween-20.
- 5) Wash 5 minutes in Tris (0.05 M, pH 7.6).
- 6) Stain 45-60 minutes with 50-100 µl primary antibody (diluted in Tris).
- 7) Wash 2 minutes in Tris.
- 8) Wash 3 minutes in Tris.
- 9) Stain 45-60 minutes with 50-100 µl secondary antibody (diluted in Tris).
- 10) Wash 2 minutes in Tris.
- 11) Wash 3 minutes in Sodium Bicarbonate buffer (0.1M, pH 8.20).
- 12) Stain 45-60 minutes with HRP-Avidin (diluted 1:5000 in Sodium Bicarbonate buffer).
- 13) Wash 2 minutes in Sodium Bicarbonate buffer.
- 14) Wash 3 minutes in Tris.
- 15) Develop for 10 minutes in DAB solution in Tris buffer.
- 16) Wash 2 minutes in Tris.
- 17) Wash 3 minutes in DPBS.
- 18) Enhance staining with 100 µl of 1% OsO₄ in DPBS for 1-5 minutes.
- 19) Wash 2 minutes in DPBS.
- 20) Post-fix slides 5 minutes in 1% Gluteraldehyde in DPBS.
- 21) Dehydrate 5 minutes in 95% alcohol.
- 22) Dehydrate 5 minutes in 100% alcohol.
- 23) Dehydrate 5 minutes in Xylene or HistoClear.
- 24) Coverslip slides using Permount.