

Mitogen Assay in Microtiter Plates
(Greg A. Perry, Ph.D.; *after Robinson et al. (1976)*)

Equipment:

Flat Bottom Microtiter Plate (Costar #3596)
Glass Fiber Filters (Brandel, or Whatman #1827-887)
PhD cell harvester (Brandel, Model 290)
Scintillation Counting Vials (Fisher #03-337-1 or Research Products Int'l #125516)

Reagents:

RF10-M Media (Sterile)
Cell Preparation (Sterile)
 - at 2×10^6 cells/ml in RF10-M media
Mitogens (pre-diluted)
 - Concanavalin A (Con-A) (Type IV; Sigma #C-2010)
 - Phytohemagglutinin-A (PHA) (Type M; Gibco #10576-015 or Sigma #L-8902)
 - Lipopolysaccharide (LPS) (Serotype 055:B5; Sigma #L-2880)
 ^3H -Thymidine (Amersham)
Scintillation Counting Fluid (PCS, Amersham #NPCS104 or Bio-Safe NA, Research Products Int'l #111198)

Method:

1. Add 100 μl of cells to each appropriate well of the microtiter plate. (2×10^5 cells per well)
2. Add 100 μl of diluted mitogen (or media) to each appropriate well containing cells.
3. Incubate at 37°C in a humidified atmosphere of 5% CO_2 in air for 52 hours.
4. Dilute thymidine to a concentration of $20 \mu\text{Ci/ml}$ in RF10-M.
5. Add 50 μl of diluted thymidine to each well (1 $\mu\text{Ci/well}$).
6. Incubate at 37°C in a humidified atmosphere of 5% CO_2 in air for another 15-18 hours.
7. Harvest cells onto glass fiber filters using a cell harvester.
8. Allow filters to air dry overnight on aluminum foil.
9. Place filter disks into liquid scintillation bags and add 2ml PCS liquid scintillation fluid.
10. Count in a liquid scintillation counter.
11. Express results as either ΔCPM or as Stimulation Index (SI) as calculated below.

Note: Each new batch of mitogen must be tittered to determine optimal concentration for stimulation.
Typical concentrations would be:

Con-A: ~ 5 $\mu\text{g/ml}$ (1-5 $\mu\text{g/ml}$)
PHA: ~ 5 $\mu\text{g/ml}$ (1-5 $\mu\text{g/ml}$)
LPS: ~ 40 $\mu\text{g/ml}$ (2-50 $\mu\text{g/ml}$)

ΔCPM : ΔCPM is calculated as ...

$$\Delta\text{CPM} = (\text{Mean}_{\text{CPM}} \text{ of mitogen stimulated wells}) - (\text{Mean}_{\text{CPM}} \text{ of media stimulated wells})$$

For example, if the mean CPM of the stimulated wells is 2350, and the mean CPM of the media stimulated wells is 124, then the $\Delta\text{CPM} = 2350 - 124 = 2226$.

Stimulation Index (SI): Stimulation index is calculated as ...

$$SI = (\text{Mean}_{\text{CPM}} \text{ of mitogen stimulated wells}) / (\text{Mean}_{\text{CPM}} \text{ of media stimulated wells})$$

If we use the same values as above, then $SI = 2350 / 124 = 18.95$

Reference: Robinson JH and JJT Owen, *Generation of T-cell function in organ culture of foetal mouse thymus: I. Mitogen responsiveness.*, *Clin. Exp. Immunol.* 23:347-354, 1976.