

Standard Technique for Mouse Cell Preparations

(Greg A. Perry, Ph.D.)

(Note: Can be used for Spleen, Thymus, Lymph Node, Peyer's Patch, and others)

Equipment:

Mouse	60mm Petri dish
Needles	Pipettes
Sterile Scissors	1cc Syringe
15ml Conical tubes (2 per tissue)	Ice bucket & ice

Reagents:

70% Ethanol
Sterile media (HBSS or culture media)
ACT (Ammonium Chloride Tris) lysis solution
Sterile media with serum (5% or more serum is fine)

Methods:

- 1) Sacrifice mouse by cervical dislocation.
- 2) Wet fur with 70% ethanol.
- 3) Clip fur and remove tissue of interest to 60mm Petri dish containing sterile media. Keep on ice.
- 4) Remove cells from the tissue by:
 - a. Teasing apart the tissue using needles.
 - b. Aspirate the tissue up and down using a 1cc syringe until most of the cells are removed from the tissue.
- 5) Decant the cells (and tissue ghost) into a 15 ml conical tube and let sit on ice for 5 minutes.
(Note: This lets the remaining tissue fragments settle out.)
- 6) Remove the supernatant into a clean 15ml conical tube.
- 7) Centrifuge for 10 minutes at 350g (approx. 1000-2000 rpm).
- 8) Discard supernatant.

Tissues containing a significant number of red blood cells (such as spleen, blood and bone marrow) must be shock-lysed prior to doing cell counts. Tissues without RBC's (such as Thymus, lymph node and Peyer's patch) do not need shock-lysing. For these tissues, proceed directly to step 13.

- 9) Resuspend cell pellet in 5ml of ACT and allow to sit on ice for 5 minutes.
- 10) Add 5ml of serum containing media (example: PBS+5%FBS).
- 11) Centrifuge for 10 minutes @ 350g (approx. 1000-2000 rpm).
- 12) Discard supernatant.
- 13) Resuspend cell pellet in known volume of media for cell counts.

<u>General guidelines on final volume:</u>	Spleen	1-2 ml (depending on size)
	Thymus	1 ml
	Lymph Node	1 ml
	Peyer's' Patch	1 ml