

Subculture of Adherent Cells
(Greg A. Perry, Ph.D.)

Equipment:

Pipettes (sterile)
Tissue culture flask (25 cm² or other size; sterile)

Reagents:

Media without serum (eg. RPMI; sterile)
Media with serum (eg. RF10; sterile)
Trypsin solution (0.25% in media without serum)
or
EDTA solution (1 mM EDTA)

Method:

- 1) Remove the old media from the flask.
- 2) Wash the flask with media without serum.
 - Add 1/2 volume of media without serum to the side of the flask opposite the cells.
 - Tip the media onto the cells to rinse them.
 - Remove and discard the rinse.
- 3) Add trypsin solution to the flask.
 - Add 3-5ml trypsin to the side of the flask opposite the cells (3-5ml / 25cm² flask).
 - Tip the trypsin onto the cells.
 - Leave the trypsin on the cells for approximately 15-30 seconds.
 - Remove and discard the trypsin.
- 4) Allow the flask to incubate until the cells "round up". Usually 5-15 minutes.
- 5) Tip the flask to allow the monolayer to slide down the surface.
- 6) Add 3 ml media (with serum) to the flask.
- 7) Resuspend the cells in the media by repeatedly gently pipetting over the surface of the flask.
- 8) Add 9 ml of fresh media with serum to 3 new 25 cm² flasks (properly labeled).
- 9) Add 1 ml of cell suspension to each of the 3 new flasks.
- 10) Cap new flasks and return to incubator.

Notes:

- This procedure assumes use of 25 cm² tissue culture flasks.
- The total amount of media in a flask should be approximately 0.2–0.5 ml / cm². More than this amount will limit gaseous diffusion through the media to the cells.

| Type | Size | Area (cm ²) | Minimum (ml) | Maximum (ml) | Standard (ml) |
|------------|---------|-------------------------|--------------|--------------|---------------|
| Flask | T-25 | 25 | 5 ml | 12.5 ml | 10 |
| | T-75 | 75 | 15 ml | 37.5 ml | 25 |
| | T-150 | 150 | 30 ml | 75 ml | 50 |
| Petri Dish | 60 mm | 21 | 4.2 | 10.5 | 8 |
| | 100 mm | 55 | 11 | 27.5 | 20 |
| Plate | 6-well | 9.5 | 1.9 | 4.75 | 3.5 |
| | 12-well | 3.8 | 0.76 | 1.9 | 1.5 |
| | 24-well | 1.9 | 0.38 | 0.95 | 0.75 |
| | 48-well | 0.8 | 0.16 | 0.4 | 0.3 |
| | 96-well | 0.32 | 0.064 | 0.16 | 0.15 |

- The activity of trypsin is negated in the presence of serum. Thus each flask should be washed with a serum free media prior to trypsinization.
- 1mM EDTA can be substituted for trypsin if the adherent cells are difficult to remove. When using EDTA, either:
 - substitute EDTA for trypsin and follow step 3 in the above protocol exactly, or ...
 - prewash flask with 3-5ml EDTA, remove, then trypsinize as in step 3 above, or ...
 - use EDTA only (5ml/25cm² flask) and leave on for 5-15 minutes, then resume at step 5.
- It should be noted that EDTA may be toxic to some cell types.